

# Ningning Wang

## List of Publications by Year in descending order

Source: <https://exaly.com/author-pdf/4750202/publications.pdf>

Version: 2024-02-01

15  
papers

472  
citations

840776

11  
h-index

996975

15  
g-index

15  
all docs

15  
docs citations

15  
times ranked

773  
citing authors

#	ARTICLE	IF	CITATIONS
1	Recent advances in chromophore-assembled upconversion nanoprobe for chemo/biosensing. <i>TrAC - Trends in Analytical Chemistry</i> , 2022, 151, 116602.	11.4	20
2	Target MicroRNA-Responsive DNA Hydrogel-Based Surface-Enhanced Raman Scattering Sensor Arrays for MicroRNA-Marked Cancer Screening. <i>Analytical Chemistry</i> , 2020, 92, 2649-2655.	6.5	78
3	Two-Photon Excitation/Red Emission, Ratiometric Fluorescent Nanoprobe for Intracellular pH Imaging. <i>Analytical Chemistry</i> , 2020, 92, 583-587.	6.5	34
4	Microsphere-based suspension array for simultaneous recognition and quantification of multiple cancer-associated miRNA via DNAzyme-Mediated signal amplification. <i>Analytica Chimica Acta</i> , 2020, 1140, 69-77.	5.4	7
5	Upconversion Nanoprobes for in Vitro and ex Vivo Measurement of Carbon Monoxide. <i>ACS Applied Materials &amp; Interfaces</i> , 2019, 11, 26684-26689.	8.0	22
6	Two-photon excitation nanoprobe for DNases activity imaging assay in hepatic ischemia reperfusion injury. <i>Sensors and Actuators B: Chemical</i> , 2019, 298, 126853.	7.8	6
7	Cyclodextrin supramolecular inclusion-enhanced pyrene excimer switching for highly selective detection of RNase H. <i>Analytica Chimica Acta</i> , 2019, 1088, 137-143.	5.4	13
8	Hybridization-activated spherical DNAzyme for cascading two-photon fluorescence emission: Applied for intracellular miRNA measurement by two-photon microscopy. <i>Sensors and Actuators B: Chemical</i> , 2019, 286, 250-257.	7.8	17
9	A spherical nucleic acid-based two-photon nanoprobe for RNase H activity assay in living cells and tissues. <i>Nanoscale</i> , 2019, 11, 8133-8137.	5.6	12
10	Oligonucleotide Cross-Linked Hydrogel for Recognition and Quantitation of MicroRNAs Based on a Portable Glucometer Readout. <i>ACS Applied Materials &amp; Interfaces</i> , 2019, 11, 7792-7799.	8.0	50
11	Polycarbonate-based core-crosslinked redox-responsive nanoparticles for targeted delivery of anticancer drug. <i>Journal of Materials Chemistry B</i> , 2018, 6, 3348-3357.	5.8	20
12	Peptide-fluorophore/AuNP conjugate-based two-photon excited fluorescent nanosensor for caspase-3 activity imaging assay in living cells and tissue. <i>MedChemComm</i> , 2017, 8, 1435-1439.	3.4	9
13	Upconversion Nanoprobes for the Ratiometric Luminescent Sensing of Nitric Oxide. <i>Journal of the American Chemical Society</i> , 2017, 139, 12354-12357.	13.7	147
14	Curing kinetics, thermal and mechanical properties of TDE-85 modified by bicyclo-benzoxazine. <i>RSC Advances</i> , 2016, 6, 90994-91001.	3.6	9
15	A novel AgNP/DNA/TPdye conjugate-based two-photon nanoprobe for GSH imaging in cell apoptosis of cancer tissue. <i>Chemical Communications</i> , 2015, 51, 16810-16812.	4.1	28