Patricia C Edwards

List of Publications by Year in descending order

Source: https://exaly.com/author-pdf/6498622/publications.pdf

Version: 2024-02-01

15 papers	1,469 citations	933447 10 h-index	996975 15 g-index
15	15	15	2370 citing authors
all docs	docs citations	times ranked	

#	Article	IF	Citations
1	Differences in SMA-like polymer architecture dictate the conformational changes exhibited by the membrane protein rhodopsin encapsulated in lipid nano-particles. Nanoscale, 2021, 13, 13519-13528.	5.6	10
2	High-mass MALDI-MS unravels ligand-mediated G protein–coupling selectivity to GPCRs. Proceedings of the National Academy of Sciences of the United States of America, 2021, 118, .	7.1	9
3	Molecular basis of β-arrestin coupling to formoterol-bound β1-adrenoceptor. Nature, 2020, 583, 862-866.	27.8	177
4	Strategic Screening and Characterization of the Visual GPCR-mini-G Protein Signaling Complex for Successful Crystallization. Journal of Visualized Experiments, 2020, , .	0.3	1
5	Molecular basis for high-affinity agonist binding in GPCRs. Science, 2019, 364, 775-778.	12.6	90
6	Crystal structure of rhodopsin in complex with a mini-G _o sheds light on the principles of G protein selectivity. Science Advances, 2018, 4, eaat7052.	10.3	65
7	Cryo-EM structure of the serotonin 5-HT1B receptor coupled to heterotrimeric Go. Nature, 2018, 558, 620-623.	27.8	183
8	Strategy for the Thermostabilization of an Agonist-Bound GPCR Coupled to a G Protein. Methods in Enzymology, 2017, 594, 243-264.	1.0	11
9	Mini-G proteins: Novel tools for studying GPCRs in their active conformation. PLoS ONE, 2017, 12, e0175642.	2.5	206
10	Complete Reversible Refolding of a G-Protein Coupled Receptor on a Solid Support. PLoS ONE, 2016, 11, e0151582.	2.5	11
11	Generation of Tetracycline-Inducible Mammalian Cell Lines by Flow Cytometry for Improved Overproduction of Membrane Proteins. Methods in Molecular Biology, 2016, 1432, 63-78.	0.9	8
12	The 2.1 \tilde{A} Resolution Structure of Cyanopindolol-Bound \hat{I}^2 1-Adrenoceptor Identifies an Intramembrane Na+ Ion that Stabilises the Ligand-Free Receptor. PLoS ONE, 2014, 9, e92727.	2.5	157
13	Structure of the large ribosomal subunit from human mitochondria. Science, 2014, 346, 718-722.	12.6	260
14	Stabilized G protein binding site in the structure of constitutively active metarhodopsin-II. Proceedings of the National Academy of Sciences of the United States of America, 2012, 109, 119-124.	7.1	226
15	Crystals of Native and Modified Bovine Rhodopsins and Their Heavy Atom Derivatives. Journal of Molecular Biology, 2004, 343, 1439-1450.	4.2	55