Sylwia Dworakowska

List of Publications by Year in descending order

Source: https://exaly.com/author-pdf/7551082/publications.pdf

Version: 2024-02-01

933447 1281871 12 599 10 11 citations h-index g-index papers 12 12 12 835 docs citations times ranked citing authors all docs

#	Article	IF	CITATIONS
1	A new way of creating cellular polyurethane materials: NIPU foams. European Polymer Journal, 2015, 66, 129-138.	5.4	125
2	Chemical and biological-based isoprene production: Green metrics. Catalysis Today, 2015, 239, 38-43.	4.4	93
3	Toward Green Atom Transfer Radical Polymerization: Current Status and Future Challenges. Advanced Science, 2022, 9, e2106076.	11.2	73
4	The role of catalysis in the synthesis of polyurethane foams based on renewable raw materials. Catalysis Today, 2014, 223, 148-156.	4.4	69
5	Size control mechanism of ZnO nanoparticles obtained in microwave solvothermal synthesis. Nanotechnology, 2018, 29, 065601.	2.6	64
6	Microwave-Assisted Synthesis of Polyols from Rapeseed Oil and Properties of Flexible Polyurethane Foams. Polymers, 2012, 4, 1462-1477.	4.5	59
7	Formulation of bio-based epoxy foams from epoxidized cardanol and vegetable oil amine. European Journal of Lipid Science and Technology, 2015, 117, 1893-1902.	1.5	48
8	Mesoporous molecular sieves containing niobium(V) as catalysts for the epoxidation of fatty acid methyl esters and rapeseed oil. Journal of Cleaner Production, 2017, 166, 901-909.	9.3	28
9	Copper mediated epoxidation of high oleic natural oils with a cumene–O2 system. Catalysis Communications, 2015, 64, 80-85.	3. 3	26
10	Crossâ€Linking of Polyesters Based on Fatty Acids. European Journal of Lipid Science and Technology, 2019, 121, 1900264.	1.5	10
11	Thiol-Ene Coupling of High Oleic Sunflower Oil towards Application in the Modification of Flexible Polyurethane Foams. Materials, 2022, 15, 628.	2.9	4

Kinetic comparison of isomeric oligo(ethylene oxide) (meth)acrylates: Aqueous polymerization of oligo(ethylene oxide) methyl ether methacrylate and methyl 2â€(oligo(ethylene oxide) methyl) Tj ETQq0 0 0 rgBT / ® xerlock 100 Tf 50 29