

A Katrin Helfer-Hungerbuehler

List of Publications by Year in descending order

Source: <https://exaly.com/author-pdf/935150/publications.pdf>

Version: 2024-02-01

11
papers

234
citations

1163117

8
h-index

1281871

11
g-index

11
all docs

11
docs citations

11
times ranked

310
citing authors

#	ARTICLE	IF	CITATIONS
1	Adeno-Associated Vector-Delivered CRISPR/SaCas9 System Reduces Feline Leukemia Virus Production In Vitro. <i>Viruses</i> , 2021, 13, 1636.	3.3	5
2	Decreased Sensitivity of the Serological Detection of Feline Immunodeficiency Virus Infection Potentially Due to Imported Genetic Variants. <i>Viruses</i> , 2019, 11, 697.	3.3	19
3	Pre-existing antibodies to candidate gene therapy vectors (adeno-associated vector serotypes) in domestic cats. <i>PLoS ONE</i> , 2019, 14, e0212811.	2.5	10
4	Putative progressive and abortive feline leukemia virus infection outcomes in captive jaguarundis (<i>Puma yagouaroundi</i>). <i>Virology Journal</i> , 2017, 14, 226.	3.4	5
5	Retroviral DNA – the silent winner: blood transfusion containing latent feline leukemia provirus causes infection and disease in naïve recipient cats. <i>Retrovirology</i> , 2015, 12, 105.	2.0	30
6	Long-term follow up of feline leukemia virus infection and characterization of viral RNA loads using molecular methods in tissues of cats with different infection outcomes. <i>Virus Research</i> , 2015, 197, 137-150.	2.2	44
7	No benefit of therapeutic vaccination in clinically healthy cats persistently infected with feline leukemia virus. <i>Vaccine</i> , 2015, 33, 1578-1585.	3.8	10
8	GAPDH Pseudogenes and the Quantification of Feline Genomic DNA Equivalents. <i>Molecular Biology International</i> , 2013, 2013, 1-7.	1.7	16
9	Quantification and molecular characterization of the feline leukemia virus A receptor. <i>Infection, Genetics and Evolution</i> , 2011, 11, 1940-1950.	2.3	6
10	Dominance of highly divergent feline leukemia virus A progeny variants in a cat with recurrent viremia and fatal lymphoma. <i>Retrovirology</i> , 2010, 7, 14.	2.0	22
11	Quantitative TaqMan® real-time PCR assays for gene expression normalisation in feline tissues. <i>BMC Molecular Biology</i> , 2009, 10, 106.	3.0	67