

Adaptation of a Mitochondrial Complex III Assay for Automated Reproducibility and Precision

Clinical Chemistry

49, 330-332

DOI: 10.1373/49.2.330

Citation Report

| # | ARTICLE | IF | CITATIONS |
|---|---|-----|-----------|
| 1 | Automated Spectrophotometric Analysis of Mitochondrial Respiratory Chain Complex Enzyme Activities in Cultured Skin Fibroblasts. <i>Clinical Chemistry</i> , 2005, 51, 2110-2116. | 1.5 | 51 |
| 2 | Molecular modeling and experimental evidence for hypericin as a substrate for mitochondrial complex III; mitochondrial photodamage as demonstrated using specific inhibitors. <i>Free Radical Biology and Medicine</i> , 2008, 45, 1581-1590. | 1.3 | 15 |
| 3 | Ataxia telangiectasia mutated influences cytochrome c oxidase activity. <i>Biochemical and Biophysical Research Communications</i> , 2011, 405, 599-603. | 1.0 | 33 |
| 4 | Enhanced Formation of Methylglyoxal-Derived Advanced Glycation End Products in Arabidopsis Under Ammonium Nutrition. <i>Frontiers in Plant Science</i> , 2018, 9, 667. | 1.7 | 32 |