## CITATION REPORT List of articles citing

Determination of porcine fibrinogen in rat and dog plasma after intraperitoneal injection of a porcine-derived fibrin glue by fluorescein-labeled assay method: Comparison with isotope-labeled assay metho

DOI: 10.1016/j.jpba.2011.08.037 Journal of Pharmaceutical and Biomedical Analysis, 2012, 57, 7-12.

Source: https://exaly.com/paper-pdf/54586389/citation-report.pdf

Version: 2024-04-25

This report has been generated based on the citations recorded by exaly.com for the above article. For the latest version of this publication list, visit the link given above.

The third column is the impact factor (IF) of the journal, and the fourth column is the number of citations of the article.

#	Paper	IF	Citations
7	Fibrin Sealant: The Only Approved Hemostat, Sealant, and Adhesive-a Laboratory and Clinical Perspective. <i>ISRN Surgery</i> , <b>2014</b> , 2014, 203943		177
6	Effects of a novel marine natural product: pyrano indolone alkaloid fibrinolytic compound on thrombolysis and hemorrhagic activities in vitro and in vivo. <i>Archives of Pharmacal Research</i> , <b>2015</b> , 38, 1530-40	6.1	12
5	In vivo study of novelly formulated porcine-derived fibrinogen as an efficient sealant. <i>Journal of Materials Science: Materials in Medicine</i> , <b>2015</b> , 26, 146	4.5	4
4	Preparation, Pharmacokinetic Profile, and Tissue Distribution Studies of a Liposome-Based Formulation of SN-38 Using an UPLC-MS/MS Method. <i>AAPS PharmSciTech</i> , <b>2016</b> , 17, 1450-1456	3.9	11
3	Development and validation of a reversed phase liquid chromatographic method with fluorescence detection for the pharmacokinetic study of a new chimeric peptide. <i>Analytical Methods</i> , <b>2016</b> , 8, 2620-2	627	1
2	Study of the degradation mechanisms of oxidized Cellulose-Derivative hemostatic agent by using model compound. <i>Journal of Biomaterials Science, Polymer Edition</i> , <b>2019</b> , 30, 580-591	3.5	3
1	A Novel Marine Pyran-Isoindolone Compound Enhances Fibrin Lysis Mediated by Single-Chain Urokinase-Type Plasminogen Activator. <b>2022</b> , 20, 495		