

Simultaneous determination of 20(*S*)-protopan
rat plasma by LC-MS/MS: application to their pharma

Biomedical Chromatography

32, e4252

DOI: 10.1002/bmc.4252

Citation Report

#	ARTICLE	IF	CITATIONS
1	LC-MS/MS determination of ginsenoside compound K and its metabolite 20(S)-protopanaxadiol in human plasma and urine: applications in a clinical study. <i>Bioanalysis</i> , 2019, 11, 365-380.	0.6	6
2	Identification of Adenylate Kinase 5 as a Protein Target of Ginsenosides in Brain Tissues Using Mass Spectrometry-Based Drug Affinity Responsive Target Stability (DARTS) and Cellular Thermal Shift Assay (CETSA) Techniques. <i>Journal of Agricultural and Food Chemistry</i> , 2022, 70, 2741-2751.	2.4	6
3	Elevated serum levels of kynurenine pathway metabolites in patients with Behçet disease. <i>Amino Acids</i> , 2022, 54, 877-887.	1.2	5